INSTRUCTIONS FOR USE en

Device for intra-abdominal pressure measurement

INTENDED USE:

Adapter for indirect intra-abdominal pressure measurement through a common bladder catheter.

DESCRIPTION - CONTENTS AND MATERIAL:

The sterile pre-assembled set consists of the following components:

- IAP adapter (barbed connector with funnel end) made of silicone
- integrated measuring line made of silicone
- large slide clamp made of ABS
- small slide clamp made of polyethylene
- female Luer lock hub made of polyamide
- attached male Luer lock protective cap made of polypropylene

DESCRIPTION - PACKAGING MATERIAL:

Sterile bag, cardboard box.

WARNING:

Do not use if individual package is damaged or has already been opened.

DISPOSAL:

Discard in accordance with public and in-house directives.

PREPARATION FOR USE:

Integration of the IAP adapter into the urinary drainage system.

⚠ ATTENTION Keep to common sterile procedures

- Open package. Connect the IAP adapter to the bladder catheter and the outlet tube of the urinary drainage unit. To do this, insert the barbed connector of the adapter into the catheter and impose the funnel end on the barbed connector of the outlet tube.
- Remove by rotation the protective cap from the proximal end of the measuring line.
- Connect the measuring line to a 3way stopcock.
- Establish connection of the 3way stopcock to the pressure transducer and the pressure module of the monitor the usual way.
- For air removal, completely prime the measuring line with saline via the 3way stopcock.
- Closing the small slide clamp avoids inadvertent emptying of the measuring line.

INSTRUCTIONS FOR PERFORMING THE MEASUREMENT:

- Close the funnel of the IAP adapter by sliding the large clamp while taking care that it is in the position intended (circumferential notch between measuring line and funnel end).
⚠ ATTENTION Connection to the urinary drainage unit is now interrupted.
- Set basic line at the height of the symphysis with the small slide clamp closed.

- For adult patients, inject 25ml sterile saline at most into the released measuring line, by means of an adequate syringe. For infant patients, apply correspondingly less volume (1 ml/kg for infants up to 20 kg). The solution applied will fill the bladder.
- Switch the 3way stopcock to occlude the path to the syringe.
⚠ ATTENTION Make sure the path towards the pressure transducer is open.
- Remove syringe.
- Read the bladder pressure value displayed on the monitor.
- Glide the large clamp at the funnel end of the IAP adapter to open.
⚠ ATTENTION Outlet of saline and urine is not possible prior to opening the clamp.

INSTRUCTIONS FOR REPEATED MEASUREMENTS:

- ⚠ ATTENTION Prior to starting a new measurement, check if the measuring line is free of air if not, repeat priming.**
- Perform further measurements as described above.

INSTRUCTIONS FOR DWELL TIME

- Dwell time of the IAP Adapter complies with the dwell time of the urinary bladder catheter.
- With replacement of the bladder catheter the IAP Adapter is recommended to be changed at the same time, however at the latest after 29 days of stay even if the bladder catheter remains beyond this time frame.

INSTRUCTIONS FOR REMOVING THE IAP ADAPTER FROM THE URINARY DRAINAGE SYSTEM:

- ⚠ WARNING Use sterile procedures.**
- Make sure the empty pressure measuring line is not connected to the 3way stopcock any more.
- **RECOMMENDATION:** Keep a new urinary drainage system ready for use.
- Disconnect the IAP adapter with the urinary drainage unit attached, by pulling the barbed connector out of the funnel end of the bladder catheter.
⚠ ATTENTION! During this procedure, it is essential to avoid direct traction onto the bladder catheter.
- Connect the prepared urinary drainage unit.

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